

Ethanol (Alcohol) (ALCO)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

23 May 2026 16:13



[Click here to view the results online](#)

TEST OVERVIEW

| | |
|--------------------------|-------------------|
| Test Name | Ethanol (Alcohol) |
| Test Code | ALCO |
| Short Description | Ethanol (Alcohol) |

OVERVIEW

| | |
|--------------------|--|
| Test Name | Ethanol (Alcohol) |
| Test Code | ALCO |
| Category | Biochemistry |
| TAT | Main Lab: 6, Hour(s) Family Site: <8hrs, <6hrs |
| Specimen(s) | 1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube |

SPECIMEN(S)

SST-Serum Separator Tube

| | |
|---------------------------------------|---|
| Specimen Type | SST-Serum Separator Tube |
| Specimen Format | Tube |
| Specimen Colour | Gold |
| Specimen Volume | 5 mL |
| Sampling Order | 2 |
| Origin | Venous blood |
| Collection time after baseline | - |
| Transport Temperature | 15-25°C |
| Accepted Other Specimens | Lithium Heparin Plasma EDTA Plasma Citratd Plasma |

| | |
|-----------------------|---|
| | Serum |
| TAT | Main Lab: 6, Hour(s) Family Site: <8hrs, <6hrs |
| Test Stability | Room Temp: 2 Week(s) 2–8°C: 6 Month(s) |

CLINICAL INFORMATION

Ethanol (Alcohol)

| | |
|---|--------------------------|
| Methodology | - |
| Specimen Type | SST-Serum Separator Tube |
| Delay before pre-treatment | 8 |
| Transport Temperature | 15-25°C |
| Transport Stability at room temp | 2 - |
| Transport Stability at 2–8°C | 6 - |
| Haemolysis interference | No |

Clinical Interest

The **blood alcohol concentration (BAC)** is a measure of the concentration of ethanol (alcohol) in the blood. It has important clinical, legal and public health applications.

High levels of alcohol in the blood can cause symptoms ranging from euphoria and disinhibition to drowsiness, confusion, respiratory depression, coma and, in extreme cases, death.

The level of alcohol in the blood determines the **seriousness of the intoxication**. For example, a BAC of over 0.3 g/L is associated with a high risk of ethyl coma, while a BAC of over 4 g/L is potentially fatal. Although BAC alone is not sufficient to diagnose chronic alcoholism, repeated high levels may indicate excessive and problematic alcohol consumption, warranting further evaluation.

In patients undergoing treatment for alcohol dependence, regular BAC testing can help monitor abstinence and identify relapses. Values and interpretation:

- Blood alcohol level of 0.1 to 0.5 g/L: Mild euphoria, disinhibition.
- Blood alcohol level of 0.5 to 1.5 g/L: Impaired coordination, slowed reflexes, impaired judgement.
- Alcohol levels of 1.5 to 3 g/L: Confusion, drowsiness, unsteady gait, vomiting.
- Alcohol levels of 3 to 4 g/L: Possible alcohol-induced coma, respiratory depression.
- Alcohol levels > 4 g/L: Risk of death from respiratory depression and cardiovascular collapse.

PATIENT INFORMATION

| | |
|--------------------------------------|---|
| Clinical Information Required | - |
| Patient Collection Notes | - |

COMMENTS & NOTES

LOINC Code

859-0, 74859-0, 5639-0

Outwork

No