

## Creatinine (24hrs Urine) (CR24)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

05 April 2026 02:50



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Creatinine (24hrs Urine)
<b>Test Code</b>	CR24
<b>Short Description</b>	Creatinine (24hrs Urine)

### OVERVIEW

<b>Test Name</b>	Creatinine (24hrs Urine)
<b>Test Code</b>	CR24
<b>Category</b>	Urine Biochemistry
<b>TAT</b>	Main Lab: 4 Hour(s), Day(s) Family Site: <4hrs
<b>Specimen(s)</b>	1 x Urine - 2000 mL 24-Hours urine container - Red - Urine 24hrs - No preservative 1 x - - 20 mL Sterile Urine container - Red - Urine Random No Preservative

### SPECIMEN(S)

#### Urine 24hrs - No preservative

<b>Specimen Type</b>	Urine 24hrs - No preservative
<b>Specimen Format</b>	24-Hours urine container
<b>Specimen Colour</b>	Red
<b>Specimen Volume</b>	2000 mL
<b>Sampling Order</b>	0
<b>Origin</b>	Urine
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	-

**TAT** Main Lab: 4 Hour(s), Day(s)  
Family Site: <4hrs

**Test Stability** Room Temp: -  
2–8°C: -

### Urine Random No Preservative

**Specimen Type** Urine Random No Preservative

**Specimen Format** Sterile Urine container

**Specimen Colour** Red

**Specimen Volume** 20 mL

**Sampling Order** 0

**Origin** -

**Collection time after baseline** -

**Transport Temperature** -

**Accepted Other Specimens** -

**TAT** Main Lab: 4 Hour(s), Day(s)  
Family Site: <4hrs

**Test Stability** Room Temp: -  
2–8°C: -

## PATIENT INFORMATION

**Clinical Information Required** 24-hour urine volume (mL)  
Collection start date  
Collection start time  
Collection end time  
Collection End date

### Patient Collection Notes

24-hour urine collection:

- On the first day, on waking, eliminate the first urine in the toilet (note the date and time).
- For the next 24 hours, collect all the urine from the day and night, including that from the next morning when you get up at the same time.
- Recap and store the container between 2 and 8 degrees C. between each micturition.
- Note the date and time of the end of micturition on the bottle.
- Be sure to close the container securely for transport.
- Bring all the urine to the laboratory as soon as possible after collection.

## COMMENTS & NOTES

**LOINC Code** 684-5, 14684-5

**Outwork**

No