

Albumin (Fluid) (FALB)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:45



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Albumin (Fluid)
Test Code	FALB
Short Description	Albumin (Fluid)

OVERVIEW

Test Name	Albumin (Fluid)
Test Code	FALB
Category	Fluid Biochemistry
TAT	Main Lab: 2 Hour(s) Family Site: <2hrs
Specimen(s)	1 x Fluid - 5 mL Steril container - Red - Fluid

SPECIMEN(S)

Fluid

Specimen Type	Fluid
Specimen Format	Steril container
Specimen Colour	Red
Specimen Volume	5 mL
Sampling Order	0
Origin	Fluid
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 2 Hour(s)

Test Stability

Family Site: <2hrs
Room Temp: 2 Month(s)
2-8°C: 5 Month(s)

PATIENT INFORMATION

Clinical Information Required	Collection_Origin
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	47-5, 1747-5
Outwork	<input type="button" value="No"/>