

## Lipase (Fluid) (FLIPA)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

15 April 2026 16:36



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Lipase (Fluid)
<b>Test Code</b>	FLIPA
<b>Short Description</b>	Lipase (Fluid)

### OVERVIEW

<b>Test Name</b>	Lipase (Fluid)
<b>Test Code</b>	FLIPA
<b>Category</b>	Fluid Biochemistry
<b>TAT</b>	Main Lab: 6 Hour(s) Family Site: <8hrs
<b>Specimen(s)</b>	1 x Fluid - 5 mL Steril container - Red - Fluid

### SPECIMEN(S)

#### Fluid

<b>Specimen Type</b>	Fluid
<b>Specimen Format</b>	Steril container
<b>Specimen Colour</b>	Red
<b>Specimen Volume</b>	5 mL
<b>Sampling Order</b>	0
<b>Origin</b>	Fluid
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	15-25°C
<b>Accepted Other Specimens</b>	-
<b>TAT</b>	Main Lab: 6 Hour(s)

**Test Stability**

Family Site: <8hrs  
Room Temp: 7 Day(s)  
2-8°C: 7 Day(s)

**PATIENT INFORMATION****Clinical Information Required**

Collection\_Origin

**Patient Collection Notes**

-

**COMMENTS & NOTES****LOINC Code**

212-4, 15212-4

**Outwork**