

Fructosamine (FRUC)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:46



[Click here to view the results online](#)

TEST OVERVIEW

| | |
|--------------------------|--------------|
| Test Name | Fructosamine |
| Test Code | FRUC |
| Short Description | FRUC |

OVERVIEW

| | |
|--------------------|--|
| Test Name | Fructosamine |
| Test Code | FRUC |
| Category | Referred tests |
| TAT | Main Lab: - - Family Site: - |
| Specimen(s) | 1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube |

SPECIMEN(S)

SST-Serum Separator Tube

| | |
|---------------------------------------|--------------------------|
| Specimen Type | SST-Serum Separator Tube |
| Specimen Format | Tube |
| Specimen Colour | Gold |
| Specimen Volume | 5 mL |
| Sampling Order | 2 |
| Origin | Venous blood |
| Collection time after baseline | - |
| Transport Temperature | - |
| Accepted Other Specimens | Serum |
| TAT | Main Lab: - - |

Test Stability

Family Site: -

Room Temp: -
2-8°C: -

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 069-8, 15069-8

Outwork Yes