

Glucose (Fasting) (GLUF)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

28 April 2026 15:03



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Glucose (Fasting)
Test Code	GLUF
Short Description	Glu (Fasting)

OVERVIEW

Test Name	Glucose (Fasting)
Test Code	GLUF
Category	Biochemistry
TAT	Main Lab: 2 Hour(s) Family Site: <2hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Grey - Sodium Fluoride Plasma

SPECIMEN(S)

Sodium Fluoride Plasma

Specimen Type	Sodium Fluoride Plasma
Specimen Format	Tube
Specimen Colour	Grey
Specimen Volume	5 mL
Sampling Order	5
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	SST-Serum Separator Tube Lithium Heparin Plasma Serum

TAT

Sodium Heparin Plasma

Main Lab: 2 Hour(s)
Family Site: <2hrs

Test Stability

Room Temp: 2 Day(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required

Fasting patient

Patient Collection Notes

Fasting for at least 8 hours, you can drink water.

COMMENTS & NOTES

LOINC Code

58-6, 1558-6

Outwork

No