

Glucose (Random) 2 Hours (GLUR2)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

07 July 2026 09:14



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Glucose (Random) 2 Hours
Test Code	GLUR2
Short Description	Glucose (Random) 2 Hours

OVERVIEW

Test Name	Glucose (Random) 2 Hours
Test Code	GLUR2
Category	Biochemistry
TAT	Main Lab: 2 Hour(s) Family Site: <3hrs, <2hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Grey - Sodium Fluoride Plasma Random 2 1 x - - 5 mL Tube - Grey - Sodium Fluorid Plama 2Hrs

SPECIMEN(S)

Sodium Fluoride Plasma Random 2

Specimen Type	Sodium Fluoride Plasma Random 2
Specimen Format	Tube
Specimen Colour	Grey
Specimen Volume	5 mL
Sampling Order	5
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-

TAT Main Lab: 2 Hour(s)
Family Site: <3hrs, <2hrs

Test Stability Room Temp: 2 Day(s)
2–8°C: 7 Day(s)

Sodium Fluorid Plasma 2Hrs

Specimen Type Sodium Fluorid Plasma 2Hrs

Specimen Format Tube

Specimen Colour Grey

Specimen Volume 5 mL

Sampling Order 5

Origin -

Collection time after baseline 2

Transport Temperature 15-25°C

Accepted Other Specimens -

TAT Main Lab: 2 Hour(s)
Family Site: <3hrs, <2hrs

Test Stability Room Temp: 2 Day(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 749-6, 14749-6

Outwork