

Glucose (Random) 2 Hours (GLUR2)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

28 April 2026 15:02



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Glucose (Random) 2 Hours
Test Code	GLUR2
Short Description	Glucose (Random) 2 Hours

OVERVIEW

Test Name	Glucose (Random) 2 Hours
Test Code	GLUR2
Category	Biochemistry
TAT	Main Lab: 2 Hour(s) Family Site: <2hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Grey - Sodium Fluoride Plasma Random 2 1 x - - 5 mL Tube - Grey - Sodium Fluorid Plama 2Hrs

SPECIMEN(S)

Sodium Fluoride Plasma Random 2

Specimen Type	Sodium Fluoride Plasma Random 2
Specimen Format	Tube
Specimen Colour	Grey
Specimen Volume	5 mL
Sampling Order	5
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-

TAT	Main Lab: 2 Hour(s) Family Site: <2hrs
Test Stability	Room Temp: 2 Day(s) 2–8°C: 7 Day(s)

Sodium Fluorid Plasma 2Hrs

Specimen Type	Sodium Fluorid Plasma 2Hrs
Specimen Format	Tube
Specimen Colour	Grey
Specimen Volume	5 mL
Sampling Order	5
Origin	-
Collection time after baseline	2
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 2 Hour(s) Family Site: <2hrs
Test Stability	Room Temp: 2 Day(s) 2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	749-6, 14749-6
Outwork	<input type="button" value="No"/>