

Indirect Anti-Globulin Test (IAT) (IAT)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:46



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Indirect Anti-Globulin Test (IAT)
Test Code	IAT
Short Description	IAT

OVERVIEW

Test Name	Indirect Anti-Globulin Test (IAT)
Test Code	IAT
Category	Immunohaematology
TAT	Main Lab: 6, Hour(s) Family Site: <6hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Plasma

SPECIMEN(S)

EDTA Plasma

Specimen Type	EDTA Plasma
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	SST-Serum Separator Tube Serum

TAT

Main Lab: 6, Hour(s)
Family Site: <6hrs

Test Stability

Room Temp: 24 Hour(s)
2–8°C: 72 Hour(s)

PATIENT INFORMATION

Clinical Information Required

Pregnancy

-

-

Patient Collection Notes

-

COMMENTS & NOTES

LOINC Code

08-2, 1008-2

Outwork

No