

Myoglobin (MYOG)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:39



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Myoglobin
Test Code	MYOG
Short Description	Myoglobin

OVERVIEW

Test Name	Myoglobin
Test Code	MYOG
Category	POC Cardiac
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs, <6hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Green - Lithium Heparin Whole Blood 1 x - - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

Lithium Heparin Whole Blood

Specimen Type	Lithium Heparin Whole Blood
Specimen Format	Tube
Specimen Colour	Green
Specimen Volume	5 mL
Sampling Order	3
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	Sodium Heparin Whole Blood Lithium Heparin Plasma

Sodium Heparin Plasma
EDTA Plasma
Serum

TAT Main Lab: 6 Hour(s)
Family Site: <8hrs, <6hrs

Test Stability Room Temp: 6 Hour(s)
2–8°C: -

SST-Serum Separator Tube

Specimen Type SST-Serum Separator Tube

Specimen Format Tube

Specimen Colour Gold

Specimen Volume 5 mL

Sampling Order 2

Origin -

Collection time after baseline -

Transport Temperature -

Accepted Other Specimens Sodium Heparin Whole Blood
Lithium Heparin Plasma
Sodium Heparin Plasma
EDTA Plasma
Serum

TAT Main Lab: 6 Hour(s)
Family Site: <8hrs, <6hrs

Test Stability Room Temp: 6 Hour(s)
2–8°C: -

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 39-3, 2639-3

Outwork