

HLA B27-PCR (HLAB27_PCR)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:34



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	HLA B27-PCR
Test Code	HLAB27_PCR
Short Description	HLA B27-PCR

OVERVIEW

Test Name	HLA B27-PCR
Test Code	HLAB27_PCR
Category	Referred tests
TAT	Main Lab: 1, 72 Day(s), Hour(s) Family Site: 1 Day
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Whole Blood

SPECIMEN(S)

EDTA Whole Blood

Specimen Type	EDTA Whole Blood
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	EDTA Whole Blood
TAT	Main Lab: 1, 72 Day(s), Hour(s)

Test Stability

Family Site: 1 Day
Room Temp: 7 Day(s)
2-8°C: 30 Day(s)

PATIENT INFORMATION**Clinical Information Required**

-

Patient Collection Notes

A separate EDTA (lavender-top) tube must be submitted for all HLA testing. If submitted with any other test requiring an EDTA tube, please submit a separate tube for this test.

COMMENTS & NOTES**LOINC Code**

043-0, 26043-0

Outwork**Yes**