

Tacrolimus (TACRO)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:37



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Tacrolimus
Test Code	TACRO
Short Description	TACRO

OVERVIEW

Test Name	Tacrolimus
Test Code	TACRO
Category	Biochemistry
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Lavender - EDTA Whole Blood

SPECIMEN(S)

EDTA Whole Blood

Specimen Type	EDTA Whole Blood
Specimen Format	Tube
Specimen Colour	Lavender
Specimen Volume	5 mL
Sampling Order	4
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-20°C
Accepted Other Specimens	-
TAT	Main Lab: 6 Hour(s)

Test Stability

Family Site: <8hrs

Room Temp: 8 Hour(s)
2–8°C: 7 Day(s)**PATIENT INFORMATION****Clinical Information Required**Date of last medication taken
Drug name
Time of last medication taken
Frequency**Patient Collection Notes**

Specimen collection should be performed before administering a new dose of Tacrolimus and always at the same time.

COMMENTS & NOTES**LOINC Code**

721-3, 32721-3

Outwork