

Toxoplasma IgG (TOXG)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:33



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Toxoplasma IgG
Test Code	TOXG
Short Description	TOXO IgG

OVERVIEW

Test Name	Toxoplasma IgG
Test Code	TOXG
Category	Serology
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	2-8°C
Accepted Other Specimens	Lithium Heparin Plasma Serum EDTA Plasma

TAT

Sodium Heparin Plasma

Main Lab: 24 Hour(s)
Family Site: 1 Day

Test Stability

Room Temp: 3 Day(s)
2–8°C: 14 Day(s)

PATIENT INFORMATION

Clinical Information Required

Pregnancy
Date of pregnancy

Patient Collection Notes

-

COMMENTS & NOTES

LOINC Code

88-4, 5388-4, 22580-5

Outwork

No