

pH (Dipstick) (UPH)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:45



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	pH (Dipstick)
Test Code	UPH
Short Description	UPH

OVERVIEW

Test Name	pH (Dipstick)
Test Code	UPH
Category	Urine Biochemistry
TAT	Main Lab: 4 Hour(s) Family Site: <4hrs
Specimen(s)	1 x Urine - 20 mL Sterile Urine container - Red - Urine Random No Preservative

SPECIMEN(S)

Urine Random No Preservative

Specimen Type	Urine Random No Preservative
Specimen Format	Sterile Urine container
Specimen Colour	Red
Specimen Volume	20 mL
Sampling Order	0
Origin	Urine
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	- Urine24H-No preservative

TAT

Main Lab: 4 Hour(s)
Family Site: <4hrs

Test Stability

Room Temp: 24 Hour(s)
2-8°C: -

PATIENT INFORMATION

Clinical Information Required

-

Patient Collection Notes

-

COMMENTS & NOTES

LOINC Code

03-2, 5803-2

Outwork

No