

Propoxyphene (Urine) (UPROP)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

04 May 2026 14:18



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Propoxyphene (Urine)
Test Code	UPROP
Short Description	Propoxyphene (Urine)

OVERVIEW

Test Name	Propoxyphene (Urine)
Test Code	UPROP
Category	Urine Drug Screen
TAT	Main Lab: 12 Hour(s) Family Site: <12hrs
Specimen(s)	1 x Urine - 20 mL Sterile Urine container - Red - Urine Random No Preservative

SPECIMEN(S)

Urine Random No Preservative

Specimen Type	Urine Random No Preservative
Specimen Format	Sterile Urine container
Specimen Colour	Red
Specimen Volume	20 mL
Sampling Order	0
Origin	Urine
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 12 Hour(s)

Test Stability

Family Site: <12hrs
Room Temp: 24 Hour(s)
2-8°C: 2 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 141-1, 19141-1

Outwork