

Anti Streptolysin O (ASOT) (ASO)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

16 June 2026 02:50



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Anti Streptolysin O (ASOT)
Test Code	ASO
Short Description	Anti-Streptolisina O (ASO)

OVERVIEW

Test Name	Anti Streptolysin O (ASOT)
Test Code	ASO
Category	Allergy/Autoimmune
TAT	Main Lab: 6, 8 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Lithium Heparin Plasma Sodium Heparin Plasma Serum

TAT

Main Lab: 6, 8 Hour(s)
Family Site: <8hrs

Test Stability

Room Temp: 1 Day(s)
2–8°C: 2 Day(s)

CLINICAL INFORMATION

Anti Streptolysin O (ASO)

Methodology

-

Specimen Type

SST-Serum Separator Tube

Delay before pre-treatment

5

Transport Temperature

15-25°C

Transport Stability at room temp

1 Day

Transport Stability at 2–8°C

2 Day

Haemolysis interference

No

PATIENT INFORMATION

Clinical Information Required

-

Patient Collection Notes

-

COMMENTS & NOTES

LOINC Code

228-4, 45228-4

Outwork

No