

Anti-centromeres IgG (CENT)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:46



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Anti-centromeres IgG
Test Code	CENT
Short Description	Anti-centromeres IgG

OVERVIEW

Test Name	Anti-centromeres IgG
Test Code	CENT
Category	Allergy/Autoimmune
TAT	Main Lab: 1 Day(s) Family Site: 1 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-20°C
Accepted Other Specimens	EDTA Plasma Sodium Heparin Plasma Serum

TAT

Main Lab: 1 Day(s)
Family Site: 1 Day(s)

Test Stability

Room Temp: 3 Hour(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 68-9, 8068-9

Outwork