

Factor IX Activity (FACT9)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:35



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Factor IX Activity
Test Code	FACT9
Short Description	FACT9

OVERVIEW

Test Name	Factor IX Activity
Test Code	FACT9
Category	Haemostasis
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Light Blue - Citrated Plasma

SPECIMEN(S)

Citrated Plasma

Specimen Type	Citrated Plasma
Specimen Format	Tube
Specimen Colour	Light Blue
Specimen Volume	5 mL
Sampling Order	1
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-20°C
Accepted Other Specimens	-
TAT	Main Lab: 24 Hour(s)

Test Stability

Family Site: 1 Day(s)
Room Temp: 3 Hour(s)
2–8°C: -

PATIENT INFORMATION**Clinical Information Required**

Hypocoagulation_investigation

Patient Collection Notes

A completely filled tube is necessary because the correct ratio of blood to citrate is critical (9:1). Mix by gentle inversion 3-4 times. Do not uncap.

COMMENTS & NOTES**LOINC Code**

87-2, 3187-2

Outwork

No