

Herpes simplex Type 2 IgG (HSV2 IgG) (HSV2G)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:45



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Herpes simplex Type 2 IgG (HSV2 IgG)
Test Code	HSV2G
Short Description	HSV2G

OVERVIEW

Test Name	Herpes simplex Type 2 IgG (HSV2 IgG)
Test Code	HSV2G
Category	Serology
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	2-8°C
Accepted Other Specimens	Serum
TAT	Main Lab: 24 Hour(s)

Test Stability

Family Site: 1 Day(s)
Room Temp: 8 Hour(s)
2–8°C: 7 Day(s)

CLINICAL INFORMATION**Herpes Simplex virus type 2 (HSV2) IgG**

Methodology	-
Specimen Type	SST-Serum Separator Tube
Delay before pre-treatment	3
Transport Temperature	2-8°C
Transport Stability at room temp	8 Hours
Transport Stability at 2–8°C	7 Day
Haemolysis interference	<input type="button" value="No"/>

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	12-9, 7912-9
Outwork	<input type="button" value="No"/>