

Treponema Pallidum Hemagglutination Assay (TPHA)

Source: Cerba Lancet Africa
View our website: cerbalancetafrica.com
Date & Time of Export: 03 April 2026 18:49



[Click here to view the results online](#)

TEST OVERVIEW

| | |
|--------------------------|---|
| Test Name | Treponema Pallidum Hemagglutination Assay |
| Test Code | TPHA |
| Short Description | TPHA |

OVERVIEW

| | |
|--------------------|--|
| Test Name | Treponema Pallidum Hemagglutination Assay |
| Test Code | TPHA |
| Category | Serology |
| TAT | Main Lab: 4 Hour(s) Family Site: <4hrs |
| Specimen(s) | 1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube |

SPECIMEN(S)

SST-Serum Separator Tube

| | |
|---------------------------------------|--|
| Specimen Type | SST-Serum Separator Tube |
| Specimen Format | Tube |
| Specimen Colour | Gold |
| Specimen Volume | 5 mL |
| Sampling Order | 2 |
| Origin | Venous blood |
| Collection time after baseline | - |
| Transport Temperature | 15-25°C |
| Accepted Other Specimens | EDTA Plasma EDTA Whole Blood Serum |

TAT

Main Lab: 4 Hour(s)
Family Site: <4hrs

Test Stability

Room Temp: 8 Day(s)
2–8°C: 3 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 032-5, 57032-5

Outwork