

Toxoplasma (IgG & IgM) (TOXO)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

15 April 2026 16:45



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Toxoplasma (IgG & IgM)
Test Code	TOXO
Short Description	TOXOP Ab

OVERVIEW

Test Name	Toxoplasma (IgG & IgM)
Test Code	TOXO
Category	Serology
TAT	Main Lab: 24 Hour(s) Family Site: 1 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	2-8°C
Accepted Other Specimens	-
TAT	Main Lab: 24 Hour(s)

Test Stability

Family Site: 1 Day(s)
Room Temp: 3 Day(s)
2-8°C: 14 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code -

Outwork