

## Glucose (Dipstick) (UGLUD)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

15 April 2026 16:46



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Glucose (Dipstick)
<b>Test Code</b>	UGLUD
<b>Short Description</b>	Glucose (Dipstick)

### OVERVIEW

<b>Test Name</b>	Glucose (Dipstick)
<b>Test Code</b>	UGLUD
<b>Category</b>	Urine Biochemistry
<b>TAT</b>	Main Lab: 4 Hour(s) Family Site: <4hrs
<b>Specimen(s)</b>	1 x Urine - 20 mL Sterile Urine container - Red - Urine Random No Preservative

### SPECIMEN(S)

#### Urine Random No Preservative

<b>Specimen Type</b>	Urine Random No Preservative
<b>Specimen Format</b>	Sterile Urine container
<b>Specimen Colour</b>	Red
<b>Specimen Volume</b>	20 mL
<b>Sampling Order</b>	0
<b>Origin</b>	Urine
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	15-25°C
<b>Accepted Other Specimens</b>	- Urine24H-No preservative

**TAT**

Main Lab: 4 Hour(s)  
Family Site: <4hrs

**Test Stability**

Room Temp: 24 Hour(s)  
2-8°C: -

### PATIENT INFORMATION

**Clinical Information Required**

-

**Patient Collection Notes**

-

### COMMENTS & NOTES

**LOINC Code**

428-4, 25428-4

**Outwork**

No