

Widal - Typhid Fever (WIDAL)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

11 May 2026 19:34



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Widal - Typhid Fever
Test Code	WIDAL
Short Description	Widal

OVERVIEW

Test Name	Widal - Typhid Fever
Test Code	WIDAL
Category	Serology
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 6 Hour(s)

Test Stability

Family Site: <8hrs
Room Temp: 24 Hour(s)
2–8°C: 1 Week(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 956-9, 53956-9

Outwork