

Creatine Kinase MB (CK-MB) (CKMB)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

07 July 2026 09:13



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Creatine Kinase MB (CK-MB)
Test Code	CKMB
Short Description	CKMB

OVERVIEW

Test Name	Creatine Kinase MB (CK-MB)
Test Code	CKMB
Category	POC Cardiac
TAT	Main Lab: 6, 2 Hour(s) Family Site: <8hrs, <6hrs, <2hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Green - Lithium Heparin Whole Blood 1 x - - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

Lithium Heparin Whole Blood

Specimen Type	Lithium Heparin Whole Blood
Specimen Format	Tube
Specimen Colour	Green
Specimen Volume	5 mL
Sampling Order	3
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	Lithium Heparin Plasma Sodium Heparin Plasma

TAT Serum
EDTA Plasma
Main Lab: 6, 2 Hour(s)
Family Site: <8hrs, <6hrs, <2hrs

Test Stability Room Temp: 4 Hour(s)
2–8°C: 8 Hour(s)

SST-Serum Separator Tube

Specimen Type SST-Serum Separator Tube

Specimen Format Tube

Specimen Colour Gold

Specimen Volume 5 mL

Sampling Order 2

Origin -

Collection time after baseline -

Transport Temperature -

Accepted Other Specimens Lithium Heparin Plasma
Sodium Heparin Plasma
Serum
EDTA Plasma

TAT Main Lab: 6, 2 Hour(s)
Family Site: <8hrs, <6hrs, <2hrs

Test Stability Room Temp: 4 Hour(s)
2–8°C: 8 Hour(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 551-5, 49551-5, 13969-1

Outwork