

## Creatine Kinase MB (CK-MB) (CKMB)

Source:

Cerba Lancet Africa

View our website:

[cerbalancetafrica.com](http://cerbalancetafrica.com)

Date & Time of Export:

28 April 2026 15:08



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Creatine Kinase MB (CK-MB)
<b>Test Code</b>	CKMB
<b>Short Description</b>	CKMB

### OVERVIEW

<b>Test Name</b>	Creatine Kinase MB (CK-MB)
<b>Test Code</b>	CKMB
<b>Category</b>	POC Cardiac
<b>TAT</b>	Main Lab: 6, 2 Hour(s) Family Site: <8hrs, <6hrs, <2hrs
<b>Specimen(s)</b>	1 x Venous blood - 5 mL Tube - Green - Lithium Heparin Whole Blood 1 x - - 5 mL Tube - Gold - SST-Serum Separator Tube

### SPECIMEN(S)

#### Lithium Heparin Whole Blood

<b>Specimen Type</b>	Lithium Heparin Whole Blood
<b>Specimen Format</b>	Tube
<b>Specimen Colour</b>	Green
<b>Specimen Volume</b>	5 mL
<b>Sampling Order</b>	3
<b>Origin</b>	Venous blood
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	Lithium Heparin Plasma Sodium Heparin Plasma

	Serum EDTA Plasma
<b>TAT</b>	Main Lab: 6, 2 Hour(s) Family Site: <8hrs, <6hrs, <2hrs
<b>Test Stability</b>	Room Temp: 4 Hour(s) 2–8°C: 8 Hour(s)

### SST-Serum Separator Tube

<b>Specimen Type</b>	SST-Serum Separator Tube
<b>Specimen Format</b>	Tube
<b>Specimen Colour</b>	Gold
<b>Specimen Volume</b>	5 mL
<b>Sampling Order</b>	2
<b>Origin</b>	-
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	Lithium Heparin Plasma Sodium Heparin Plasma Serum EDTA Plasma
<b>TAT</b>	Main Lab: 6, 2 Hour(s) Family Site: <8hrs, <6hrs, <2hrs
<b>Test Stability</b>	Room Temp: 4 Hour(s) 2–8°C: 8 Hour(s)

### PATIENT INFORMATION

<b>Clinical Information Required</b>	-
<b>Patient Collection Notes</b>	-

### COMMENTS & NOTES

<b>LOINC Code</b>	551-5, 49551-5, 13969-1
<b>Outwork</b>	<input type="button" value="No"/>