

Copper (COPPER)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:35



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Copper
Test Code	COPPER
Short Description	Copp

OVERVIEW

Test Name	Copper
Test Code	COPPER
Category	Biochemistry
TAT	Main Lab: 7 Day(s) Family Site: 7 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Red - Serum

SPECIMEN(S)

Serum

Specimen Type	Serum
Specimen Format	Tube
Specimen Colour	Red
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 7 Day(s)

Test Stability

Family Site: 7 Day(s)
Room Temp: 5 Day(s)
2–8°C: 10 Day(s)

PATIENT INFORMATION**Clinical Information Required**

-

Patient Collection Notes

Draw blood into a trace metal tube (BD Royal Blue cap).

COMMENTS & NOTES**LOINC Code**

665-4, 14665-4

Outwork**Yes**