

Osmolality Urine (OSMOLALITY_URINE)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:33



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Osmolality Urine
Test Code	OSMOLALITY_URINE
Short Description	Uosmal

OVERVIEW

Test Name	Osmolality Urine
Test Code	OSMOLALITY_URINE
Category	Referred tests
TAT	Main Lab: 1 Day(s) Family Site: 1 Day
Specimen(s)	1 x Urine - 20 mL Sterile Urine container - Red - Urine Random No Preservative

SPECIMEN(S)

Urine Random No Preservative

Specimen Type	Urine Random No Preservative
Specimen Format	Sterile Urine container
Specimen Colour	Red
Specimen Volume	20 mL
Sampling Order	0
Origin	Urine
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Urine Urine24H-No preservative

TAT

Main Lab: 1 Day(s)
Family Site: 1 Day

Test Stability

Room Temp: 4 Day(s)
2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 95-5, 2695-5

Outwork Yes