

Malignant Melanoma S-100B (S-100B)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:35



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Malignant Melanoma S-100B
Test Code	S-100B
Short Description	Malignant Melanoma S-100B

OVERVIEW

Test Name	Malignant Melanoma S-100B
Test Code	S-100B
Category	Referred tests
TAT	Main Lab: 24 Hour(s) Family Site: -
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-20°C
Accepted Other Specimens	Serum
TAT	Main Lab: 24 Hour(s)

Test Stability

Family Site: -
Room Temp: 12 Hour(s)
2-8°C: 48 Hour(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 275-3, 47275-3

Outwork Yes