

RAST F14 Soya Bean (RAST F14 SOYA)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:31



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	RAST F14 Soya Bean
Test Code	RAST F14 SOYA
Short Description	RAST F14 Soya Bean

OVERVIEW

Test Name	RAST F14 Soya Bean
Test Code	RAST F14 SOYA
Category	Referred tests
TAT	Main Lab: 72 Hour(s) Family Site: 3 Days
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 72 Hour(s)

Test Stability

Family Site: 3 Days
Room Temp: 14 Day(s)
2-8°C: 14 Day(s)

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 568-9, 15568-9

Outwork Yes