

Aspergillus Galactomannan Antigen (ASPERG_GALAC)

Source: Cerba Lancet Africa
View our website: cerbalancetafrica.com
Date & Time of Export: 27 June 2026 10:40



[Click here to view the results online](#)

TEST OVERVIEW

| | |
|--------------------------|-----------------------------------|
| Test Name | Aspergillus Galactomannan Antigen |
| Test Code | ASPERG_GALAC |
| Short Description | Asperg Galac |

OVERVIEW

| | |
|--------------------|--|
| Test Name | Aspergillus Galactomannan Antigen |
| Test Code | ASPERG_GALAC |
| Category | Referred tests |
| TAT | Main Lab: 7 Day(s) Family Site: - |
| Specimen(s) | 1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube |

SPECIMEN(S)

SST-Serum Separator Tube

| | |
|---------------------------------------|--------------------------|
| Specimen Type | SST-Serum Separator Tube |
| Specimen Format | Tube |
| Specimen Colour | Gold |
| Specimen Volume | 5 mL |
| Sampling Order | 2 |
| Origin | Venous blood |
| Collection time after baseline | - |
| Transport Temperature | 15-25°C |
| Accepted Other Specimens | Serum |
| TAT | Main Lab: 7 Day(s) |

Test Stability

Family Site: -
Room Temp: 24 Hour(s)
2-8°C: -

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 383-9, 35383-9

Outwork Yes