

Anti-cardiolipin antibody (aCL IgG, A &M) (ANTI_CARDIOLIPIN)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:35



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Anti-cardiolipin antibody (aCL IgG, A &M)
Test Code	ANTI_CARDIOLIPIN
Short Description	Anti-cardiolipin antibody (aCL IgG, A &M)

OVERVIEW

Test Name	Anti-cardiolipin antibody (aCL IgG, A &M)
Test Code	ANTI_CARDIOLIPIN
Category	Allergy/Autoimmune
TAT	Main Lab: - - Family Site: -
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum

TAT

Main Lab: - -
Family Site: -

Test Stability

Room Temp: 5 Day(s)
2-8°C: 14 Day(s)

CLINICAL INFORMATION

Anti-cardiolipin antibody (aCL IgG, A &M)

Methodology -

Specimen Type SST-Serum Separator Tube

Delay before pre-treatment 6

Transport Temperature 15-25°C

Transport Stability at room temp 5 Day

Transport Stability at 2-8°C 14 Day

Haemolysis interference

PATIENT INFORMATION

Clinical Information Required -

Patient Collection Notes -

COMMENTS & NOTES

LOINC Code 62-2, 8062-2

Outwork