

SARS Covid-19 Antigen (COVAG)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

11 June 2026 01:48



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	SARS Covid-19 Antigen
Test Code	COVAG
Short Description	COVAG

OVERVIEW

Test Name	SARS Covid-19 Antigen
Test Code	COVAG
Category	Molecular biology
TAT	Main Lab: 1 Hour(s) Family Site: <1hr
Specimen(s)	1 x Nose - - mL - - - - Nasopharyngeal swab

SPECIMEN(S)

Nasopharyngeal swab

Specimen Type	Nasopharyngeal swab
Specimen Format	-
Specimen Colour	-
Specimen Volume	- mL
Sampling Order	-
Origin	Nose
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	-
TAT	Main Lab: 1 Hour(s)

Test Stability

Family Site: <1hr
Room Temp: 2 Hour(s)
2-8°C: -

CLINICAL INFORMATION**SARS Covid-19 Antigen**

Methodology	Manual
Specimen Type	Nasopharyngeal swab
Delay before pre-treatment	30
Transport Temperature	15-25°C
Transport Stability at room temp	2 Hours
Transport Stability at 2-8°C	- -
Haemolysis interference	<input type="button" value="No"/>

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	558-4, 94558-4
Outwork	<input type="button" value="No"/>