

Anti-Factor Xa(Clexane) (ANTI -FACTOR Xa)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:24



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Anti-Factor Xa(Clexane)
Test Code	ANTI -FACTOR Xa
Short Description	Anti-Factor Xa

OVERVIEW

Test Name	Anti-Factor Xa(Clexane)
Test Code	ANTI -FACTOR Xa
Category	Haemostasis
TAT	Main Lab: 7 Day(s) Family Site: 7 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Light Blue - Citrated Plasma

SPECIMEN(S)

Citrated Plasma

Specimen Type	Citrated Plasma
Specimen Format	Tube
Specimen Colour	Light Blue
Specimen Volume	5 mL
Sampling Order	1
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	-
TAT	Main Lab: 7 Day(s)

Test Stability

Family Site: 7 Day(s)

Room Temp: -
2-8°C: -**CLINICAL INFORMATION****Anti-Factor Xa(Clexane)**

Methodology	-
Specimen Type	Citrated Plasma
Delay before pre-treatment	-
Transport Temperature	-
Transport Stability at room temp	--
Transport Stability at 2-8°C	--
Haemolysis interference	<input type="button" value="No"/>

PATIENT INFORMATION

Clinical Information Required -
Patient Collection Notes

The citrated tube must be collected 3 hours after the Clexane Injection was given.

COMMENTS & NOTES

LOINC Code -
Outwork