LDL Cholesterol (Direct) (LDLD)

Source:

View our website:

Date & Time of Export:

Click here to view the results online

Cerba Lancet Africa cerbalancetafrica.com

05 November 2025 19:10



TEST OVERVIEW

Test Name

LDL Cholesterol (Direct)

Test Code

LDLD

Short Description

LDL Cholesterol (Direct)

OVERVIEW

Test Name

LDL Cholesterol (Direct)

Test Code

LDLD

Category

Biochemistry

TAT

Main Lab: 6 Hour(s) Family Site: <8hrs

Specimen(s)

1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type SST-Serum Separator Tube

 Specimen Format
 Tube

 Specimen Colour
 Gold

 Specimen Volume
 5 mL

 Sampling Order
 2

Origin Venous blood

Collection time after baseline

Transport Temperature 15-25°C

Accepted Other Specimens Serum

Lithium Heparin Plasma

EDTA Plasma

Sodium Heparin Plasma

TAT Main Lab: 6 Hour(s)

Family Site: <8hrs

Room Temp: 8 Hour(s) 2-8°C: 5 Day(s) **Test Stability**

PATIENT INFORMATION

Clinical Information Required

Fasting patient Monitoring as part of treatment

Patient Collection Notes

In the absence of known hypertriglyceridaemia, a lipid profile may be taken in non-fasting patients. If hypertriglyceridaemia is discovered, a control sample should be prescribed and taken after a 12-hour fast. In all cases, the patient may drink water.

COMMENTS & NOTES

LOINC Code

419-0, 69419-0

Outwork