

Calcium (24 Hrs Urine) (CA24)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:33



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Calcium (24 Hrs Urine)
Test Code	CA24
Short Description	Calcium (24 Hrs Urine)

OVERVIEW

Test Name	Calcium (24 Hrs Urine)
Test Code	CA24
Category	Urine Biochemistry
TAT	Main Lab: 4 Hour(s) Family Site: <5hrs, <4hrs
Specimen(s)	1 x Urine - 2000 mL 24-Hours urine container - Red - Urine 24hrs - HCL

SPECIMEN(S)

Urine 24hrs - HCL

Specimen Type	Urine 24hrs - HCL
Specimen Format	24-Hours urine container
Specimen Colour	Red
Specimen Volume	2000 mL
Sampling Order	0
Origin	Urine
Collection time after baseline	-
Transport Temperature	-
Accepted Other Specimens	-
TAT	Main Lab: 4 Hour(s)

Test Stability

Family Site: <5hrs, <4hrs

Room Temp: -
2-8°C: -**PATIENT INFORMATION****Clinical Information Required**24-hour urine volume (mL)
Collection start date
Collection start time
Collection end time
Collection End date**Patient Collection Notes**

24-hour urine collection:

- On the first day, on waking, eliminate the first urine in the toilet (note the date and time).
- For the next 24 hours, collect all the urine from the day and night, including that from the next morning when you get up at the same time.
- Recap and store the container between 2 and 8 degrees C. between each micturition.
- Note the date and time of the end of micturition on the bottle.
- Be sure to close the container securely for transport.
- Bring all the urine to the laboratory as soon as possible after collection.

COMMENTS & NOTES**LOINC Code**

362-5, 25362-5

Outwork

No