

Chromogranin A (CHROMOA)

Source:

Cerba Lancet Africa

View our website:

cerbalancetafrica.com

Date & Time of Export:

27 June 2026 10:34



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Chromogranin A
Test Code	CHROMOA
Short Description	Chromogranin A

OVERVIEW

Test Name	Chromogranin A
Test Code	CHROMOA
Category	Referred tests
TAT	Main Lab: 7 Day(s) Family Site: 7 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	-20°C
Accepted Other Specimens	Serum
TAT	Main Lab: 7 Day(s)

Test Stability

Family Site: 7 Day(s)
Room Temp: 7 Day(s)
2–8°C: - Day(s)

PATIENT INFORMATION**Clinical Information Required**

Date of last medication taken
-

Patient Collection Notes

To be taken at least 15 days before starting treatment with proton pump inhibitors (PPIs).

COMMENTS & NOTES**LOINC Code**

169-7, 30169-7

Outwork

Yes