

## Contactin-associated protein 2 IgG Ab (CASPR2-Ab)

Source: Cerba Lancet Africa  
View our website: [cerbalancetafrica.com](http://cerbalancetafrica.com)  
Date & Time of Export: 15 April 2026 16:37



[Click here to view the results online](#)

### TEST OVERVIEW

<b>Test Name</b>	Contactin-associated protein 2 IgG Ab
<b>Test Code</b>	CASPR2-Ab
<b>Short Description</b>	CASPR2-Ab

### OVERVIEW

<b>Test Name</b>	Contactin-associated protein 2 IgG Ab
<b>Test Code</b>	CASPR2-Ab
<b>Category</b>	Referred tests
<b>TAT</b>	Main Lab: - - Family Site: -
<b>Specimen(s)</b>	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

### SPECIMEN(S)

#### SST-Serum Separator Tube

<b>Specimen Type</b>	SST-Serum Separator Tube
<b>Specimen Format</b>	Tube
<b>Specimen Colour</b>	Gold
<b>Specimen Volume</b>	5 mL
<b>Sampling Order</b>	2
<b>Origin</b>	Venous blood
<b>Collection time after baseline</b>	-
<b>Transport Temperature</b>	-
<b>Accepted Other Specimens</b>	Serum -

**TAT**

Main Lab: - -  
Family Site: -

**Test Stability**

Room Temp: -  
2-8°C: -

### PATIENT INFORMATION

**Clinical Information Required** -

**Patient Collection Notes** -

### COMMENTS & NOTES

**LOINC Code** 285-4, 94285-4

**Outwork**  Yes