

Phospholipase A2 Receptor (PLA2R) Antibody (PLA2R)

Source: Cerba Lancet Africa
View our website: cerbalancetafrica.com
Date & Time of Export: 22 May 2026 16:18



[Click here to view the results online](#)

TEST OVERVIEW

Test Name	Phospholipase A2 Receptor (PLA2R) Antibody
Test Code	PLA2R
Short Description	Phospholipase A2 Receptor (PLA2R) Antibody

OVERVIEW

Test Name	Phospholipase A2 Receptor (PLA2R) Antibody
Test Code	PLA2R
Category	Allergy/Autoimmune
TAT	Main Lab: 7 Day(s) Family Site: 7 Day(s)
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST-Serum Separator Tube

SPECIMEN(S)

SST-Serum Separator Tube

Specimen Type	SST-Serum Separator Tube
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 7 Day(s)

Test Stability

Family Site: 7 Day(s)
Room Temp: 7 Day(s)
2-8°C: 14 Day(s)

CLINICAL INFORMATION**Phospholipase A2 Receptor (PLA2R) Antibody**

Methodology	-
Specimen Type	SST-Serum Separator Tube
Delay before pre-treatment	4
Transport Temperature	15-25°C
Transport Stability at room temp	7 Day
Transport Stability at 2-8°C	14 Day
Haemolysis interference	<input type="button" value="No"/>

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	737-9, 73737-9
Outwork	<input type="button" value="Yes"/>