# Glucose (Random) 2 Hours (GLUR2)

Source: Cerba Lancet Africa

View our website: cerbalancetafrica.com

Date & Time of Export: 05 November 2026 04:45

Click here to view the results online

## TEST OVERVIEW

| | |
|---|---|
| **Test Name** | Glucose (Random) 2 Hours |
| **Test Code** | GLUR2 |
| **Short Description** | Glucose (Random) 2 Hours |

## OVERVIEW

| | |
|---|---|
| **Test Name** | Glucose (Random) 2 Hours |
| **Test Code** | GLUR2 |
| **Category** | Biochemistry |
| **TAT** | Main Lab: 2 Hour(s)<br>Family Site: <3hrs, <2hrs |
| **Specimen(s)** | 1 x Venous blood - 5 mL Tube - Grey - Sodium Fluoride Plasma Random 2<br>1 x - - 5 mL Tube - Grey - Sodium Fluorid Plama 2Hrs |

## SPECIMEN(S)

### Sodium Fluoride Plasma Random 2

| | |
|---|---|
| **Specimen Type** | Sodium Fluoride Plasma Random 2 |
| **Specimen Format** | Tube |
| **Specimen Colour** | Grey |
| **Specimen Volume** | 5 mL |
| **Sampling Order** | 5 |
| **Origin** | Venous blood |
| **Collection time after baseline** | - |
| **Transport Temperature** | 15-25°C |
| **Accepted Other Specimens** | - |

| | |
|---|---|
| **TAT** | Main Lab: 2 Hour(s)<br>Family Site: <3hrs, <2hrs |
| **Test Stability** | Room Temp: 2 Day(s)<br>2–8°C: 7 Day(s) |

### Sodium Fluorid Plama 2Hrs

| | |
|---|---|
| **Specimen Type** | Sodium Fluorid Plama 2Hrs |
| **Specimen Format** | Tube |
| **Specimen Colour** | Grey |
| **Specimen Volume** | 5 mL |
| **Sampling Order** | 5 |
| **Origin** | - |
| **Collection time after baseline** | 2 |
| **Transport Temperature** | 15-25°C |
| **Accepted Other Specimens** | - |
| **TAT** | Main Lab: 2 Hour(s)<br>Family Site: <3hrs, <2hrs |
| **Test Stability** | Room Temp: 2 Day(s)<br>2–8°C: 7 Day(s) |

## PATIENT INFORMATION

| | |
|---|---|
| **Clinical Information Required** | - |
| **Patient Collection Notes** | - |

## COMMENTS & NOTES

| | |
|---|---|
| **LOINC Code** | 749-6, 14749-6 |
| **Outwork** | No |