Lipase (LIPA)

Source: Cerba Lancet Africa

View our website: cerbalancetafrica.com

Date & Time of Export: 05 November 2026 04:43

Click here to view the results online

TEST OVERVIEW

Test Name	Lipase
Test Code	LIPA
Short Description	Lipase

OVERVIEW

Test Name	Lipase
Test Code	LIPA
Category	Biochemistry
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs, <6hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST

SPECIMEN(S)

SST

Specimen Type	SST
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Lithium Heparin Plasma Sodium Heparin Plasma Serum

TAT	Main Lab: 6 Hour(s) Family Site: <8hrs, <6hrs
Test Stability	Room Temp: 6 Day(s) 2–8°C: 7 Day(s)

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	40-3, 3040-3
Outwork	No