Widal - Typhid Fever (WIDAL)

Source: Cerba Lancet Africa

View our website: cerbalancetafrica.com

Date & Time of Export: 07 November 2026 18:46

Click here to view the results online

TEST OVERVIEW

Test Name	Widal - Typhid Fever
Test Code	WIDAL
Short Description	Widal

OVERVIEW

Test Name	Widal - Typhid Fever
Test Code	WIDAL
Category	Serology
TAT	Main Lab: 6 Hour(s) Family Site: <8hrs
Specimen(s)	1 x Venous blood - 5 mL Tube - Gold - SST

SPECIMEN(S)

SST

Specimen Type	SST
Specimen Format	Tube
Specimen Colour	Gold
Specimen Volume	5 mL
Sampling Order	2
Origin	Venous blood
Collection time after baseline	-
Transport Temperature	15-25°C
Accepted Other Specimens	Serum
TAT	Main Lab: 6 Hour(s)

Test Stability	Family Site: <8hrs Room Temp: 24 Hour(s) 2–8°C: 1 Week(s)

PATIENT INFORMATION

Clinical Information Required	-
Patient Collection Notes	-

COMMENTS & NOTES

LOINC Code	956-9, 53956-9
Outwork	No